

# Switch Content for Popular Scanning Software

IntelliTools software	Control + Enter (one-switch) Mouse Click (one-switch) Space, 0, 1, 2, 3 (multiple switch)
Speaking Dynamically Pro (Mayer Johnson)	Mouse Click (one switch) Any defined keystroke (one-switch or multiple switch)
Laureate Software	Mouse Click (one switch) 1,2 or Enter (one Switch)
Edmark Software	Space (one-Switch) Mouse Click (one-switch)
Inclusive Software	Mouse Click (one-switch) 1,2 Enter or Space (one-switch) 1 and 2 (two-switch) Enter and Space (two-switch)
Crick Software (Clicker 4/ClozePro)  Clicker 5	Left Arrow (one-switch) Up and Left Arrow (two-switch) F7 (one-switch) F7 and F8 (two-switch)
Widgit Software (Writing with Symbols)	Left Arrow (one-switch) Up and Left Arrow (two-switch)
SoftTouch Software	Mouse Click
Attainment	Space (one-switch) 1,2,3 (multiple switch)
Don Johnston	Click (one-switch) 1 (one-switch)
Linda Burkhart	(whatever Ipics Studio wants – F12+Control+Enter) Control + Right Arrow/ Control + Down Arrow
Madentec	Click (one-switch)
Applied Human Factors (Reach Interface)	Click (one-switch) Left Click + Right Click (two-switch)
RJ Cooper	Click (one-switch) Left-Click + Right Click (two-switch)
Wivik	F11, F12